

## Antibodies Product Datasheet

<b>Product Name:</b>	GTF2H1 mouse mAb
<b>Catalog Number:</b>	N0827
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IP,ICC
<b>Dilutions:</b>	wb 1:300 icc 1:50
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human GTF2H1 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	gtf2h1
<b>Protein Name:</b>	
<b>Other Name:</b>	Basic transcription factor 2 62 kDa subunit;Basic transcription factor 62 kDa subunit; BTF 2;BTF2;BTF2 p62;General transcription factor IIH polypeptide 1;General transcription factor IIH subunit 1;General transcription factor IIH, polypeptide 1, 62kDa;gtf
<b>Human Swiss Prot No.:</b>	P32780
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of GTF2H1 and does not cross-react with related proteins.
<b>Storage and Stability:</b>	-20°C/1 year