

Antibodies Product Datasheet

Product Name:	Anti-Human Glypican-3 (ABT068) Mouse Monoclonal Antibody
Catalog Number:	V0363
Clone ID:	ABT068
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	GPC3
Gene Name:	P51654
Other Name:	DGSX, Glypican proteoglycan 3, Glypican-3 [Precursor], Gpc3, GPC3_HUMAN, GTR2 2, GTR2-2, Heparan sulphate proteoglycan, Intestinal protein OCI 5, Intestinal protein OCI-5, MXR7, OCI 5, OCI-5, OCI5, SDYS, Secreted glypican-3, SGB, SGBS, SGBS1
Isotype	IgG1, Kappa
Localization:	Hepatocellular carcinoma
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Glypican-3 protein.
Storage and Stability:	-20°C/1 year