

Antibodies Product Datasheet

Product Name:	IL-31R α Antibody
Catalog Number:	C20838
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the Internal region of human IL-31R α .
Concentration:	1 mg/ml
Gene Name:	IL31RA
Protein Name:	Interleukin-31 receptor subunit alpha
Other Name:	IL31RA; CRL3; GPL; Interleukin-31 receptor subunit alpha; IL-31 receptor subunit alpha; IL-31R subunit alpha; IL-31R-alpha; IL-31RA; Cytokine receptor-like 3; GLM-R; hGLM-R; Gp130-like monocyte receptor; Gp130-like receptor; ZcytoR17
Human Swiss Prot No.:	Q8NI17
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-31R α Polyclonal Antibody detects endogenous levels of IL-31R α protein.
Storage and Stability:	-20°C/1 year