

Antibodies Product Datasheet

Product Name:	HRT2 Antibody
Catalog Number:	C11818
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HEY2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	HEY2
Protein Name:	Hairy/enhancer-of-split related with YRPW motif protein 2
Other Name:	HEY2; BHLHB32; CHF1; GRL; HERP; HERP1; HRT2; Hairy/enhancer-of-split related with YRPW motif protein 2; Cardiovascular helix-loop-helix factor 1; hCHF1; Class B basic helix-loop-helix protein 32; bHLHb32; HES-related repressor protein 2; Ha
Human Swiss Prot No.:	Q9UBP5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HRT2 Polyclonal Antibody detects endogenous levels of HRT2 protein.
Storage and Stability:	-20°C/1 year