

Antibodies Product Datasheet

Product Name:	Blimp-1 Monoclonal Antibody
Catalog Number:	N0062
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Blimp-1 expressed in E. Coli.
Concentration:	
Gene Name:	PRDM1
Protein Name:	PR domain zinc finger protein 1
Other Name:	PRDM1; BLIMP1; PR domain zinc finger protein 1; BLIMP-1; Beta-interferon gene positive regulatory domain I-binding factor; PR domain-containing protein 1; Positive regulatory domain I-binding factor 1; PRDI-BF1; PRDI-binding factor 1
Human Swiss Prot No.:	O75626
Formulation:	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Specificity:	Blimp-1 Monoclonal Antibody detects endogenous levels of Blimp-1 protein.
Storage and Stability:	-20°C/1 year