

Antibodies Product Datasheet

Product Name:	Ephrin-A5 Antibody
Catalog Number:	L0256
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Ephrin-A5 . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	EFNA5
Protein Name:	Ephrin-A5
Other Name:	EFNA5; EPLG7; LERK7; Ephrin-A5; AL-1; EPH-related receptor tyrosine kinase ligand 7; LERK-7
Human Swiss Prot No.:	P52803
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ephrin-A5 Polyclonal Antibody detects endogenous levels of Ephrin-A5 protein.
Storage and Stability:	-20°C/1 year