

Antibodies Product Datasheet

Product Name:	SNAPC 19 Antibody
Catalog Number:	C10201
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SNAPC5. AA range:10-59
Concentration:	1 mg/ml
Gene Name:	SNAPC5
Protein Name:	snRNA-activating protein complex subunit 5
Other Name:	SNAPC5; SNAP19; snRNA-activating protein complex subunit 5; SNAPc subunit 5; Small nuclear RNA-activating complex polypeptide 5; snRNA-activating protein complex 19 kDa subunit; SNAPc 19 kDa subunit
Human Swiss Prot No.:	O75971
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SNAPC 19 Polyclonal Antibody detects endogenous levels of SNAPC 19 protein.
Storage and Stability:	-20°C/1 year