

Antibodies Product Datasheet

Product Name:	Prdx1 (Phospho-Tyr194) Antibody
Catalog Number:	A9659
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Prdx1 (Tyr194)
Concentration:	1 mg/ml
Gene Name:	PRDX1 PAGA PAGB TDPX2
Protein Name:	Prdx1 (Tyr194) , Peroxiredoxin-1
Other Name:	Peroxiredoxin-1 (EC 1.11.1.15) (Natural killer cell-enhancing factor A) (NKEF-A) (Proliferation-associated gene protein) (PAG) (Thioredoxin peroxidase 2) (Thioredoxin-dependent peroxide reductase 2)
Human Swiss Prot No.:	Q06830
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Prdx1 (phospho-Tyr194)
Storage and Stability:	-20°C/1 year