

Antibodies Product Datasheet

Product Name:	JAK2 Mouse mAb(2F6)
Catalog Number:	N1490
Reactivity:	Human, Mouse,Rat
Applications:	IHC
Dilutions:	IHC 1:50-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human JAK2
Concentration:	1 mg/ml
Gene Name:	JAK2
Protein Name:	JAK2
Other Name:	Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2)
Human Swiss Prot No.:	O60674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.63% sodium azide.
Specificity:	This antibody detects endogenous levels of JAK2 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year