

Antibodies Product Datasheet

Product Name:	ADAM 12 (Cleaved-Glu208) Antibody
Catalog Number:	L0413
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ADAM 12 (Cleaved-Glu208)
Concentration:	1 mg/ml
Gene Name:	ADAM12 MLTN UNQ346/PRO545
Protein Name:	ADAM 12 (Cleaved-Glu208)
Other Name:	Disintegrin and metalloproteinase domain-containing protein 12 (ADAM 12) (EC 3.4.24.-) (Meltrin-alpha)
Human Swiss Prot No.:	O43184
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse ADAM 12 (Cleaved-Glu208, protein was cleaved amino acid sequence between 207-208)
Storage and Stability:	-20°C/1 year