

Antibodies Product Datasheet

Product Name:	BCL-10 (Phospho-Ser138) Antibody
Catalog Number:	A8249
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BCL-10 (Phospho-Ser138)
Concentration:	1 mg/ml
Gene Name:	BCL10 CIPER CLAP
Protein Name:	BCL-10 (Phospho-Ser138)
Other Name:	B-cell lymphoma/leukemia 10 (B-cell CLL/lymphoma 10) (Bcl-10) (CARD-containing molecule enhancing NF-kappa-B) (CARD-like apoptotic protein) (hCLAP) (CED-3/ICH-1 prodomain homologous E10-like regulator) (CIPER) (Cellular homolog of vCARMEN) (cCARMEN) (Cell
Human Swiss Prot No.:	O95999
Formulation:	Liquid in PBS containing 50% glycerol, and 0.165% sodium azide.
Specificity:	This antibody detects endogenous levels of BCL-10 (Phospho-Ser138) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year