
Antibodies Product Datasheet

Product Name:	Eif4b Antibody
Catalog Number:	C42180
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Eif4b
Concentration:	1 mg/ml
Gene Name:	EIF4B
Protein Name:	Eif4b
Other Name:	Eukaryotic translation initiation factor 4B (eIF-4B); YM33117
Human Swiss Prot No.:	P23588
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.77% sodium azide.
Specificity:	This antibody detects endogenous levels of Eif4b at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year