

Antibodies Product Datasheet

Product Name:	EDRF1 Antibody
Catalog Number:	C49876
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human EDRF1
Concentration:	1 mg/ml
Gene Name:	EDRF1 C10orf137
Protein Name:	Erythroid differentiation-related factor 1
Other Name:	
Human Swiss Prot No.:	Q3B7T1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of EDRF1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year