

## Antibodies Product Datasheet

<b>Product Name:</b>	Napsin A Antibody
<b>Catalog Number:</b>	C30873
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	IHC-p 1:50-200, ELISA 1:10000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthetic peptide from human protein at AA range: 191-240
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	NAPSA NAP1 NAPA
<b>Protein Name:</b>	Napsin-A (EC 3.4.23.-) (Aspartyl protease 4) (ASP4) (Asp 4) (Napsin-1) (TA01/TA02)
<b>Other Name:</b>	Napsin-A (EC 3.4.23.-) (Aspartyl protease 4) (ASP4) (Asp 4) (Napsin-1) (TA01/TA02)
<b>Human Swiss Prot No.:</b>	O96009
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	The antibody detects endogenous Napsin A
<b>Storage and Stability:</b>	-20°C/1 year