

Antibodies Product Datasheet

Product Name:	MMP-10 Antibody
Catalog Number:	C30917
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human MMP-10. at AA range: 351-400
Concentration:	1 mg/ml
Gene Name:	MMP10
Protein Name:	MMP-10
Other Name:	Stromelysin-2 (SL-2) (EC 3.4.24.22) (Matrix metalloproteinase-10) (MMP-10) (Transin-2)
Human Swiss Prot No.:	P09238
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MMP-10 Polyclonal Antibody detects endogenous levels of MMP-10
Storage and Stability:	-20°C/1 year