

Antibodies Product Datasheet

Product Name:	Tie2 (Phospho-Ser1119) Antibody
Catalog Number:	A8193
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Tie2 (Ser1119)
Concentration:	1 mg/ml
Gene Name:	TEK TIE2 VMCM VMCM1
Protein Name:	Tie2 (Ser1119)
Other Name:	Angiopoietin-1 receptor (EC 2.7.10.1) (Endothelial tyrosine kinase) (Tunica interna endothelial cell kinase) (Tyrosine kinase with Ig and EGF homology domains-2) (Tyrosine-protein kinase receptor TEK) (Tyrosine-protein kinase receptor TIE-2) (hTIE2) (p140)
Human Swiss Prot No.:	Q02763
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Tie2 (phospho-Ser1119)
Storage and Stability:	-20°C/1 year