

Antibodies Product Datasheet

Product Name:	PPGB (20k, Cleaved-Met327) Antibody
Catalog Number:	L0471
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human PPGB (20k, Cleaved-Met327)
Concentration:	1 mg/ml
Gene Name:	CTSA PPGB
Protein Name:	PPGB (20k, Cleaved-Met327)
Other Name:	Lysosomal protective protein (EC 3.4.16.5) (Carboxypeptidase C) (Carboxypeptidase L) (Cathepsin A) (Protective protein cathepsin A) (PPCA) (Protective protein for beta-galactosidase) [Cleaved into: Lysosomal protective protein 32 kDa chain; Lysosomal prot
Human Swiss Prot No.:	P10619
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse PPGB (20k, Cleaved-Met327, protein was cleaved amino acid sequence between 326-327)
Storage and Stability:	-20°C/1 year