

Antibodies Product Datasheet

Product Name:	Ghrelin Receptor Antibody
Catalog Number:	C40192
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Ghrelin Receptor
Concentration:	1 mg/ml
Gene Name:	GHSR
Protein Name:	Growth hormone secretagogue receptor type 1 (GHS-R) (GH-releasing peptide receptor) (GHRP) (Ghrelin receptor)
Other Name:	Growth hormone secretagogue receptor type 1 (GHS-R) (GH-releasing peptide receptor) (GHRP) (Ghrelin receptor); YM3557
Human Swiss Prot No.:	Q92847
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ghrelin Receptor protein(A232) detects endogenous levels of Ghrelin Receptor
Storage and Stability:	-20°C/1 year