

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | SLC38A5 Antibody  |
| <b>Catalog Number:</b>        | C50646  |
| <b>Reactivity:</b>            | Human, Mouse, Rat   |
| <b>Applications:</b>          | WB  |
| <b>Dilutions:</b>             | WB 1:500-2000   |
| <b>Modification:</b>          | Unmodified/Total  |
| <b>Source:</b>                | Polyclonal, Rabbit,IgG  |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human SLC38A5  |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | SLC38A5 JM24 SN2 SNAT5 PP7194   |
| <b>Protein Name:</b>          | Sodium-coupled neutral amino acid transporter 5 (Solute carrier family 38 member 5) (System N transporter 2)          |
| <b>Other Name:</b>            |   |
| <b>Human Swiss Prot No.:</b>  | Q8WUX1  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | This antibody detects endogenous levels of SLC38A5 at Human, Mouse,Rat  |
| <b>Storage and Stability:</b> | Store at -20° C/1 year  |