

Antibodies Product Datasheet

Product Name:	FBXW5 Antibody
Catalog Number:	C43756
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FBXW5
Concentration:	1 mg/ml
Gene Name:	FBXW5 FBW5 PP3971
Protein Name:	FBXW5
Other Name:	F-box/WD repeat-containing protein 5 (F-box and WD-40 domain-containing protein 5)
Human Swiss Prot No.:	Q969U6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.132% sodium azide.
Specificity:	This antibody detects endogenous levels of FBXW5 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year