

Antibodies Product Datasheet

Product Name:	Oct-3/4 Antibody
Catalog Number:	C30870
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Oct-3/4 at AA range: 201-250
Concentration:	1 mg/ml
Gene Name:	POU5F1
Protein Name:	Oct-3/4
Other Name:	POU domain, class 5, transcription factor 1 (Octamer-binding protein 3) (Oct-3) (Octamer-binding protein 4) (Oct-4) (Octamer-binding transcription factor 3) (OTF-3)
Human Swiss Prot No.:	Q01860
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Oct-3/4 Polyclonal Antibody detects endogenous levels of Oct-3/4
Storage and Stability:	-20°C/1 year