Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Bcl-x (Phospho-Thr115) Antibody

Catalog Number: A8243

Reactivity: Human, Mouse, Rat

Applications: IHC-p, IF(paraffin section), ELISA

Dilutions: Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in

other applications.

Modification: Phospho Specific

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human BCL-XL around the phosphorylation site of Thr115. AA range:81-

130

Concentration: 1 mg/ml

Gene Name: BCL2L1

Protein Name: Bcl-2-like protein 1

Other Name: BCL2L1; BCL2L; BCLX; Bcl-2-like protein 1; Bcl2-L-1; Apoptosis

regulator Bcl-X

Human Swiss Prot No.: Q07817

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Phospho-Bcl-x (T115) Polyclonal Antibody detects endogenous levels of

Bcl-x protein only when phosphorylated at T115.

Storage and Stability: -20°C/1 year