

Antibodies Product Datasheet

Product Name:	Brk Antibody
Catalog Number:	B0783
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Breast Tumor Kinase. AA range:402-451
Concentration:	1 mg/ml
Gene Name:	PTK6
Protein Name:	Protein-tyrosine kinase 6
Other Name:	PTK6; BRK; Protein-tyrosine kinase 6; Breast tumor kinase; Tyrosine-protein kinase BRK
Human Swiss Prot No.:	Q13882
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Brk Polyclonal Antibody detects endogenous levels of Brk protein.
Storage and Stability:	-20°C/1 year