

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: SIRP-α1 Antibody

Catalog Number: C0322

Reactivity: Human, Mouse, Rat

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000. Not yet tested in other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Sirp alpha1. AA range:451-500

Concentration: 1 mg/ml

Gene Name: SIRPA

Protein Name: Tyrosine-protein phosphatase non-receptor type substrate 1

Other Name: SIRPA; BIT; MFR; MYD1; PTPNS1; SHPS1; SIRP; Tyrosine-protein

phosphatase non-receptor type substrate 1; SHP substrate 1; SHPS-1; Brain Ig-like molecule with tyrosine-based activation motifs; Bit; CD172

antigen-like family member A; Inhibito

Human Swiss Prot No.: P78324

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: SIRP-α1 Polyclonal Antibody detects endogenous levels of SIRP-α1

protein.

Storage and Stability: -20°C/1 year