

Antibodies Product Datasheet

Product Name:	SLC25A23 Antibody
Catalog Number:	C49880
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC25A23
Concentration:	1 mg/ml
Gene Name:	SLC25A23 APC2 MCSC2 SCAMC3
Protein Name:	Calcium-binding mitochondrial carrier protein SCaMC-3 (Mitochondrial ATP-Mg/Pi carrier protein 2) (Mitochondrial Ca(2+)-dependent solute carrier protein 2) (Small calcium-binding mitochondrial carrier protein 3) (Solute carrier family 25 member 23)
Other Name:	
Human Swiss Prot No.:	Q9BV35
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC25A23 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year