

Antibodies Product Datasheet

Product Name:	FZD9 Antibody
Catalog Number:	R12-2788
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	FZD9 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FZD9.
Concentration:	1 mg/ml
Gene Name:	FZD9
Other Name:	Frizzled-9; Fz-9; hFz9; FzE6; CD349 antigen; FZD9; FZD3;
Human Swiss Prot No.:	O00144
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	FZD9 Antibodyantibody detects endogenous levels of FZD9.
Storage and Stability:	Stable at -20°C for at least 1 year.