

Antibodies Product Datasheet

Product Name:	Rb (Acetyl-K873/K874) Antibody
Catalog Number:	D12183
Reactivity:	Human:K873+K874, Mouse:K866+K867, Rat:K865+K866
Applications:	WB, ELISA
Dilutions:	wb dilution 1:1000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized Acetyl peptide derived from human Rb. at AA range: K873/872
Concentration:	1 mg/ml
Gene Name:	RB1
Protein Name:	Rb
Other Name:	Retinoblastoma-associated protein (p105-Rb) (pRb) (Rb) (pp110)
Human Swiss Prot No.:	P06400
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Rb (Acetyl-K873/K874), It doesn't react with total protein.
Storage and Stability:	-20°C/1 year