
Antibodies Product Datasheet

Product Name:	NUB1 Antibody
Catalog Number:	C0360
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NYREN18. AA range:566-615
Concentration:	1 mg/ml
Gene Name:	NUB1
Protein Name:	NEDD8 ultimate buster 1
Other Name:	NUB1; NYREN18; NEDD8 ultimate buster 1; Negative regulator of ubiquitin-like proteins 1; Renal carcinoma antigen NY-REN-18
Human Swiss Prot No.:	Q9Y5A7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NUB1 Polyclonal Antibody detects endogenous levels of NUB1 protein.
Storage and Stability:	-20°C/1 year