

Antibodies Product Datasheet

Product Name:	SMAD7 Antibody
Catalog Number:	C42033
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	SMAD7 MADH7 MADH8
Protein Name:	Mothers against decapentaplegic homolog 7 (MAD homolog 7) (Mothers against DPP homolog 7) (Mothers against decapentaplegic homolog 8) (MAD homolog 8) (Mothers against DPP homolog 8) (SMAD family member 7) (SMAD 7) (Smad7) (hSMAD7)
Other Name:	
Human Swiss Prot No.:	O15105
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	SMAD7 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year