

Antibodies Product Datasheet

Product Name:	KIF3B Antibody
Catalog Number:	C42257
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human KIF3B
Concentration:	1 mg/ml
Gene Name:	KIF3B KIAA0359
Protein Name:	Kinesin-like protein KIF3B (HH0048) (Microtubule plus end-directed kinesin motor 3B)
Other Name:	
Human Swiss Prot No.:	O15066
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human KIF3B
Storage and Stability:	-20°C/1 year