
Antibodies Product Datasheet

Product Name:	HNF4alpha Antibody
Catalog Number:	R12-2884
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	HNF4alpha Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HNF4alpha.
Concentration:	1 mg/ml
Gene Name:	HNF4A
Other Name:	HNF-4; HNF-4-alpha; HNF4; HNF4A; Hepatocyte nuclear factor 4-alpha; NR2A1; TCF14; Transcription factor 14; Transcription factor HNF-4
Human Swiss Prot No.:	P41235
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	HNF4alpha Antibodyantibody detects endogenous levels of HNF4alpha.
Storage and Stability:	Stable at -20°C for at least 1 year.