

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Rabbit Anti-Guinea Pig IgG (H&L) - Texas Red
<b>Catalog Number:</b>	SA0544
<b>Reactivity:</b>	Guinea Pig
<b>Conjugate:</b>	Texas Red
<b>Applications:</b>	ELISA, IHC, IF, FCM
<b>Dilutions:</b>	IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration)
<b>Modification:</b>	Not Modification Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
<b>Concentration:</b>	
<b>Other Name:</b>	
<b>Formulation:</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
<b>Storage and Stability:</b>	-20°C/1 year