

Antibodies Product Datasheet

Product Name:	HABP2 (50k heavy chain, Cleaved-Arg313) Antibody
Catalog Number:	L0400
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HABP2 (50k heavy chain, Cleaved-Arg313)
Concentration:	1 mg/ml
Gene Name:	HABP2 HGFAL PHBP
Protein Name:	HABP2 (50k heavy chain, Cleaved-Arg313)
Other Name:	Hyaluronan-binding protein 2 (EC 3.4.21.-) (Factor VII-activating protease) (Factor seven-activating protease) (FSAP) (Hepatocyte growth factor activator-like protein) (Plasma hyaluronan-binding protein) [Cleaved into: Hyaluronan-binding protein 2 50 kDa
Human Swiss Prot No.:	Q14520
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat HABP2 (50k heavy chain, Cleaved-Arg313, protein was cleaved amino acid sequence between 313-314)
Storage and Stability:	-20°C/1 year