

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Vimentin (ABT019) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0475
<b>Clone ID:</b>	ABT019
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC, WB
<b>Application Dilutions</b>	IHC 1:100-200, WB 1:500-1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	VIM
<b>Gene Name:</b>	P08670
<b>Other Name:</b>	CTRCT30, Epididymis luminal protein 113, FLJ36605, HEL113, VIM, VIME_HUMAN, Vimentin
<b>Isotype</b>	IgG2b, Kappa
<b>Localization:</b>	Tonsil/ Liver
<b>Positive Tissue Control:</b>	Cytoplasmic
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Vimentin protein. In western blotting of HEK293 cell lysate, the antibody can label a 54 kDa band corresponding to Vimentin.
<b>Storage and Stability:</b>	-20°C/1 year