

Antibodies Product Datasheet

Product Name:	Tom20 Antibody
Catalog Number:	C16678
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TOMM20. AA range:15-64
Concentration:	1 mg/ml
Gene Name:	TOMM20
Protein Name:	Mitochondrial import receptor subunit TOM20 homolog
Other Name:	TOMM20; KIAA0016; Mitochondrial import receptor subunit TOM20 homolog; Mitochondrial 20 kDa outer membrane protein; Outer mitochondrial membrane receptor Tom20
Human Swiss Prot No.:	Q15388
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Tom20 Polyclonal Antibody detects endogenous levels of Tom20 protein.
Storage and Stability:	-20°C/1 year