

Antibodies Product Datasheet

Product Name:	BSSP4 Antibody
Catalog Number:	C42792
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	PRSS22 BSSP4 PRSS26 SP001LA UNQ302/PRO343
Protein Name:	Brain-specific serine protease 4 (BSSP-4) (EC 3.4.21.-) (Serine protease 22) (Serine protease 26) (Trypsin epsilon)
Other Name:	
Human Swiss Prot No.:	Q9GZN4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	BSSP4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year