

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | PRC1 (Phospho-Thr481) Antibody |
| Catalog Number: | P12-1050 |
| Reactivity: | Human |
| Applications: | WB, ELISA |
| Dilutions: | PRC1 (Phospho-Thr481) Antibody |
| Modification: | Phospho Specific |
| Source: | Rabbit Polyclonal |
| Purification: | The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide. |
| Immunogen: | The antiserum was produced against synthesized phospho-peptide derived from human PRC1 around the phosphorylation site of Thr481. |
| Concentration: | 1 mg/ml |
| Gene Name: | PRC1 |
| Other Name: | Protein regulator of cytokinesis 1 |
| Human Swiss Prot No.: | O43663 |
| Formulation: | Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Specificity: | PRC1 (Phospho-Thr481) antibody detects endogenous levels of PRC1 only when phosphorylated at Thr481. |
| Storage and Stability: | Stable at -20°C for at least 1 year. |