
Antibodies Product Datasheet

Product Name:	Cyclin D1 Monoclonal Antibody
Catalog Number:	N0175
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cyclin D1 expressed in E. Coli.
Concentration:	
Gene Name:	CCND1
Protein Name:	G1/S-specific cyclin-D1
Other Name:	CCND1; BCL1; PRAD1; G1/S-specific cyclin-D1; B-cell lymphoma 1 protein; BCL-1; BCL-1 oncogene; PRAD1 oncogene
Human Swiss Prot No.:	P24385
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cyclin D1 Monoclonal Antibody detects endogenous levels of Cyclin D1 protein.
Storage and Stability:	-20°C/1 year