

Antibodies Product Datasheet

Product Name:	CaMKII β / γ / δ (Phospho Thr287) Monoclonal Antibody(4H2)
Catalog Number:	N1554
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of CaMKII β / γ / δ (Phospho Thr287)
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	
Other Name:	CAMK2B; CAM2; CAMK2; CAMKB; Calcium/calmodulin-dependent protein kinase type II subunit beta; CaM kinase II subunit beta; CaMK-II subunit beta; CAMK2G; CAMK; CAMK-II; CAMKG; Calcium/calmodulin-dependent protein kinase type II subunit gamma;
Human Swiss Prot No.:	Q13554/Q13555/Q13557
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CaMKII β / γ / δ (Phospho Thr287) protein detects endogenous levels of CaMKII β / γ / δ (Phospho Thr287)
Storage and Stability:	-20°C/1 year