

Antibodies Product Datasheet

Product Name:	C1orf187 Antibody
Catalog Number:	C50585
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1orf187
Concentration:	1 mg/ml
Gene Name:	DRAXIN C1orf187 PSEC0258 UNQ3119/PRO10268
Protein Name:	Draxin (Dorsal inhibitory axon guidance protein) (Dorsal repulsive axon guidance protein) (Neucrin)
Other Name:	
Human Swiss Prot No.:	Q8NBI3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of C1orf187 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year