
Antibodies Product Datasheet

Product Name:	PARP Monoclonal Antibody(M3)
Catalog Number:	N1396
Reactivity:	Human,Chicken(testedbyourcustomer)
Applications:	WB
Dilutions:	WB: 1:1000-3000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of PARP
Concentration:	
Gene Name:	PARP1
Protein Name:	Poly [ADP-ribose] polymerase 1
Other Name:	Poly [ADP-ribose] polymerase 1 (PARP-1) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 1) (ARTD1) (NAD(+) ADP-ribosyltransferase 1) (ADPRT 1) (Poly[ADP-ribose] synthase 1)
Human Swiss Prot No.:	P09874
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous PARP protein.
Storage and Stability:	-20°C/1 year