

Antibodies Product Datasheet

Product Name:	IL-1 α Antibody
Catalog Number:	R12-2203
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-1 α . AA range:8-57
Concentration:	1 mg/ml
Gene Name:	IL1A
Protein Name:	Interleukin-1 alpha
Other Name:	IL1A; IL1F1; Interleukin-1 alpha; IL-1 alpha; Hematopoietin-1
Human Swiss Prot No.:	P01583
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1 α Polyclonal Antibody detects endogenous levels of IL-1 α protein.
Storage and Stability:	-20°C/1 year