

Antibodies Product Datasheet

Product Name:	IκB β Mouse mAb(3E4)
Catalog Number:	N1472
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human IκB β
Concentration:	1 mg/ml
Gene Name:	NFKBIB IKBB TRIP9
Protein Name:	IκB β
Other Name:	NF-kappa-B inhibitor beta (NF-kappa-BIB) (I-kappa-B-beta) (IκB-B) (IκB-beta) (IκappaBbeta) (Thyroid receptor-interacting protein 9) (TR-interacting protein 9) (TRIP-9)
Human Swiss Prot No.:	Q15653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.45% sodium azide.
Specificity:	This antibody detects endogenous levels of IκB β at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year