

## Antibodies Product Datasheet

|                        |                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>   | His-tag mouse mAb                                                                                                  |
| <b>Catalog Number:</b> | N1032                                                                                                              |
| <b>Reactivity:</b>     |                                                                                                                    |
| <b>Applications:</b>   | WB                                                                                                                 |
| <b>Dilutions:</b>      | wb dilution 1:1000                                                                                                 |
| <b>Modification:</b>   | Unmodified/ Total                                                                                                  |
| <b>Source:</b>         | Monoclonal Mouse                                                                                                   |
| <b>Purification:</b>   | The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>      | Synthetic peptide corresponding to fragment (HHHHHH) His Tag                                                       |
| <b>Concentration:</b>  | 1 mg/ml                                                                                                            |
| <b>Gene Name:</b>      | His-Tag                                                                                                            |
| <b>Protein Name:</b>   | His Tag                                                                                                            |
| <b>Other Name:</b>     | 6 His epitope tag;Hexa His tag;HHHHHH epitope tag;HHHHHH tag;His tag;Polyhistidine Tag.                            |

### Human Swiss Prot No.:

**Formulation:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

**Specificity:** This antibody specifically recognizes the 6X His-Tag at C-terminal or N-terminal of recombinant proteins, which can be used to detect His-Tag fusion protein of various expression vector.

**Storage and Stability:** -20°C/1 year