

Antibodies Product Datasheet

Product Name:	EphA3 Antibody
Catalog Number:	C10508
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EPHA3. AA range:831-880
Concentration:	1 mg/ml
Gene Name:	EPHA3
Protein Name:	Ephrin type-A receptor 3
Other Name:	EPHA3; ETK; ETK1; HEK; TYRO4; Ephrin type-A receptor 3; EPH-like kinase 4; EK4; hEK4; HEK; Human embryo kinase; Tyrosine-protein kinase TYRO4; Tyrosine-protein kinase receptor ETK1; Eph-like tyrosine kinase 1
Human Swiss Prot No.:	P29320
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EphA3 Polyclonal Antibody detects endogenous levels of EphA3 protein.
Storage and Stability:	-20°C/1 year