

Antibodies Product Datasheet

Product Name:	RRP1 Antibody
Catalog Number:	C49871
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human RRP1
Concentration:	1 mg/ml
Gene Name:	RRP1 D21S2056E NNP1 NOP52 RRP1A
Protein Name:	Ribosomal RNA processing protein 1 homolog A (Novel nuclear protein 1) (NNP-1) (Nucleolar protein Nop52) (RRP1-like protein)
Other Name:	
Human Swiss Prot No.:	P56182
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of RRP1 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year