

Antibodies Product Datasheet

Product Name:	EFNA3 (Cleaved-Gly214) Antibody
Catalog Number:	L0399
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EFNA3 (Cleaved-Gly214)
Concentration:	1 mg/ml
Gene Name:	EFNA3 EFL2 EPLG3 LERK3
Protein Name:	EFNA3 (Cleaved-Gly214)
Other Name:	Ephrin-A3 (EFL-2) (EHK1 ligand) (EHK1-L) (EPH-related receptor tyrosine kinase ligand 3) (LERK-3)
Human Swiss Prot No.:	P52797
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EFNA3 (Cleaved-Gly214, protein was cleaved amino acid sequence between 214-215)
Storage and Stability:	-20°C/1 year