
Antibodies Product Datasheet

Product Name:	Anti-Mammaglobin [ABT-MGB]
Catalog Number:	V0091
Clone ID:	ABT-MGB
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	Q13296
Gene Name:	SCGB2A2
Other Name:	Mammaglobin-A; Mammaglobin-1; Secretoglobin family 2A member 2; SCGB2A2; MGB1; UGB2
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Mammaglobin [ABT-MGB] detects endogenous levels of total Mammaglobin protein.
Storage and Stability:	-20°C/1 year