
Antibodies Product Datasheet

Product Name:	Vimentin (Phospho-Ser39) Antibody
Catalog Number:	A8277
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human Vimentin (Ser39)
Concentration:	1 mg/ml
Gene Name:	VIM
Protein Name:	Vimentin (Ser39)
Other Name:	Vimentin
Human Swiss Prot No.:	P08670
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human Mouse Rat Vimentin (phospho-Ser39)
Storage and Stability:	-20°C/1 year