
Antibodies Product Datasheet

Product Name:	GFP-Tag Antibody
Catalog Number:	C40014
Reactivity:	Species independent
Applications:	WB, IF/ICC,
Dilutions:	WB: 1:5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of GFP-Tag
Concentration:	
Gene Name:	
Protein Name:	GFP Tag
Other Name:	YM3206
Human Swiss Prot No.:	
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects GFP and GFP tag fusion proteins.
Storage and Stability:	-20°C/1 year